



Study Agent Transport Container

Using a suitable transport container is important to ensure that the study agent is not released into the environment or result in an exposure during transportation.

IMPORTANT FEATURES

Hard-sided. Biohazard Sticker Affixed. Shatter-proof. Non-porous. Sealable Lid.

EXAMPLES

Specimen Sample Transport Box, Tupperware™ Container, Small Cooler, Tackle Box



TIPS

- Better to have a horizontal than vertical orientation
- Place an absorbent material on the bottom of the container

UNACCEPTABLE

Tall Containers, Glass Containers, Soft-sided Lunch bag/cooler, Styrofoam

